

Product Datasheet

MPT64 Antibody - FITC Conjugated (orb2226464)

| | |
|----------------------------|---|
| Description | MPT64 Antibody - FITC Conjugated |
| Species/Host | Rabbit |
| Conjugation | FITC |
| Tested Applications | WB |
| Immunogen | Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA) |
| Preservatives | Liquid |
| Form/Appearance | Liquid |
| Storage | -20°C or -80°C |
| Note | For research use only |
| Isotype | IgG |
| Clonality | Polyclonal |
| Purity | Caprylic Acid Ammonium Sulfate Precipitation purified |
| Uniprot ID | P9WIN9 |
| Expiration Date | 12 months from date of receipt. |