

Product Datasheet

MPT64 Antibody - HRP Conjugated [Out of stock] (orb2226465)

Description	MPT64 Antibody - HRP Conjugated [Out of stock]
Species/Host	Rabbit
Conjugation	HRP
Tested Applications	ELISA
Immunogen	Recombinant Mycobacterium tuberculosis Immunogenic protein MPT64 protein (24-228AA)
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Caprylic Acid Ammonium Sulfate Precipitation purified
Uniprot ID	P9WIN9
Expiration Date	12 months from date of receipt.