

Product Datasheet

YWHAH Antibody - Biotin Conjugated [Out of stock] (orb2227208)

Description	YWHAH Antibody - Biotin Conjugated [Out of stock]
Species/Host	Rabbit
Conjugation	Biotin
Tested Applications	ELISA, IHC
Immunogen	Recombinant Human 14-3-3 protein eta protein (4-246AA)
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	Caprylic Acid Ammonium Sulfate Precipitation purified
Uniprot ID	Q04917
NCBI	NP_003396.1
Expiration Date	12 months from date of receipt.