



Product Datasheet

CD101 antibody (orb222746)


Description	Mouse monoclonal antibody to CD101
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, IHC-Fr, IP, WB
Immunogen	Human thymic clone B12
Target	CD101
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Store at 2-8°C. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q93033
Entrez	9398
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml.
Expiration Date	12 months from date of receipt.



Flow cytometry multicolor surface staini...



Separation of human CD101 positive CD16 ...



Flow cytometry surface staining pattern ...