


Product Datasheet

DDIT4L antibody (orb222759)

Description	Mouse monoclonal antibody to DDIT4L
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC, WB
Immunogen	N-terminal recombinant fragment of human DDIT4L (amino acids 2-98)
Target	DDIT4L
Preservatives	Phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	1 mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified by protein-A affinity chromatography.
Uniprot ID	Q96D03
Entrez	115265
Dilution Range	Flow cytometry: Recommended dilution: 1-4 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Western blotting analysis of DDIT4L expr...