


Product Datasheet

DDIT4L antibody (PE) (orb222760)

Description	Mouse monoclonal antibody against DDIT4L conjugated to PE
Reactivity	Human
Conjugation	PE
Tested Applications	FC
Immunogen	N-terminal recombinant fragment of human DDIT4L (amino acids 2-98)
Target	DDIT4L
Preservatives	Stabilizing phosphate buffered saline (PBS), pH 7.4, 15 mM sodium azide
Concentration	0.1 mg/ml
Storage	Store at 2-8°C. Protect from prolonged exposure to light. Do not freeze.
Note	For research use only
Application notes	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Isotype	Mouse IgG1
Clonality	Monoclonal
Purity	Purified antibody is conjugated with R-phycoerythrin (PE) under optimum conditions. Unconjugated antibody and free fluorochrome are removed by size-exclusion chromatography.
Uniprot ID	Q96D03
Entrez	115265
Dilution Range	Flow cytometry: Recommended dilution: 1-5 µg/ml. Intracellular staining.
Expiration Date	12 months from date of receipt.



Separation of DDIT4L transfected HEK-293...