

Product Datasheet

Mouse anti- Sulfametoxydiazine (orb2227690)

Description	Mouse anti- Sulfametoxydiazine
Species/Host	Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Sulfametoxydiazine-BSA conjugate
Preservatives	Liquid
Form/Appearance	Liquid
Storage	-20°C or -80°C
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	Caprylic Acid Ammonium Sulfate Precipitation purified
Expiration Date	12 months from date of receipt.