


## Product Datasheet

### Human Transthyretin protein (orb222833)

<b>Description</b>	Human Transthyretin protein
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	SDS-PAGE, WB
<b>Target</b>	Transthyretin
<b>Form/Appearance</b>	0.14mol/L NaCl, 0.027mol/L KCl, 0.01mol/L Na <sub>2</sub> PO <sub>4</sub> , 0.018mol/L KH <sub>2</sub> PO <sub>4</sub> , pH 7.4
<b>Concentration</b>	1.0 mg/mL
<b>Storage</b>	Store at 4°C for short terms and -20°C for long terms storage
<b>Note</b>	For research use only
<b>Protein Sequence</b>	MASHRLLLLC LAGLVFVSEA GPTGTGESKC PLMVKVLDVAV RGSPAINVAV HVFRKAADD T WEPFASGKTS ESGELHGLTT EEEFVEGIYK VEIDTKSYWK ALGISPFHEH AEVVFTANDS GPRRYTIAAL LSPYSYSTTA VVTNPKE
<b>Purity</b>	> 85% as determined by SDS-PAGE
<b>MW</b>	21 kDa
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human TTR protein