



Product Datasheet

EGFR (Phospho-Thr678) antibody (orb222835)

Description	Rabbit polyclonal antibody to EGFR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P00533
Dilution Range	WB: 1:300-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:40000
Expiration Date	12 months from date of receipt.



Immunofluorescence analysis of HuvEc cel...



Western blot analysis of cos7 cells usin...