


Product Datasheet

JNK1/2/3 (Phospho-Thr183+Tyr185) antibody (orb222866)

| | |
|----------------------------|--|
| Description | Rabbit polyclonal antibody to JNK1/2/3 |
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | Unconjugated |
| Tested Applications | WB |
| Immunogen | Internal |
| Concentration | 1mg/ml |
| Storage | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles. |
| Note | For research use only |
| Clonality | Polyclonal |
| Uniprot ID | P45984, P53779, P45983 |
| Dilution Range | WB: 1:500-1:3000, ELISA: 1:10000 |
| Expiration Date | 12 months from date of receipt. |



Western blot analysis of 293 cells using...