



Product Datasheet

AKT (Phospho-Thr308) antibody (orb222952)

Description	Rabbit polyclonal antibody to AKT (Phospho-Thr308)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P31749, P31751
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:10000
Expiration Date	12 months from date of receipt.



Western blot analysis of 293 cells treated with insulin.



Immunohistochemical analysis of lung carcinoma tissue.