
Product Datasheet

**Human EGFR Antibody : PE (Phospho-Tyr1068)
(orb2229599)**

Description

Human EGFR Antibody : PE

| | |
|----------------------------|---|
| Species/Host | Rabbit |
| Reactivity | Human, Mouse, Rat |
| Conjugation | PE |
| Tested Applications | FC |
| Immunogen | A synthetic phospho-peptide corresponding to residues surrounding Tyr1068 of human phospho EGFR. |
| Preservatives | Liquid. 1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin. |
| Form/Appearance | Liquid. 1X PBS, 0.09% NaN ₃ , 0.2% rAlbumin. |
| Storage | 2-8C |
| Note | For research use only |
| Application notes | Application Info: For flow cytometric staining, the suggested use of this reagent is 5 uL per million cells or 5 uL per 100 uL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application. |
| Isotype | Rabbit IgG1k |
| Clonality | Monoclonal |
| Uniprot ID | P00533 |
| NCBI | NP_005219.2 |
| Expiration Date | 12 months from date of receipt. |