

---

## Product Datasheet

**Human EGFR Antibody : PE (Phospho-Tyr1068)  
(orb2229600)**

**Description**

Human EGFR Antibody : PE

<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	PE
<b>Tested Applications</b>	FC
<b>Immunogen</b>	A synthetic phospho-peptide corresponding to residues surrounding Tyr1068 of human phospho EGFR.
<b>Preservatives</b>	Liquid. 1X PBS, 0.09% NaN <sub>3</sub> , 0.2% rAlbumin.
<b>Form/Appearance</b>	Liquid. 1X PBS, 0.09% NaN <sub>3</sub> , 0.2% rAlbumin.
<b>Storage</b>	2-8C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: For flow cytometric staining, the suggested use of this reagent is 5 uL per million cells or 5 uL per 100 uL of staining volume. It is recommended that the reagent be titrated for optimal performance for each application.
<b>Isotype</b>	Rabbit IgG1k
<b>Clonality</b>	Monoclonal
<b>Uniprot ID</b>	<a href="#">P00533</a>
<b>NCBI</b>	<a href="#">NP_005219.2</a>
<b>Expiration Date</b>	12 months from date of receipt.