



Product Datasheet

FKHR (Phospho-Ser256) antibody (orb223004)

Description	Rabbit polyclonal antibody to FKHR
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q12778
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using FKHR/FOXO1A (pSer256) antibody.



Immunohistochemical analysis of human brain tissue using FKHR/FOXO1A (pSer256) antibody.