



Product Datasheet

IRS-1 (Phospho-Ser312) antibody (orb223032)


Description	Rabbit polyclonal antibody to IRS-1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, WB
Immunogen	Internal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P35568
Dilution Range	WB: 1:500-1:3000, IF/ICC: 1:100-1:500, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HepG2 cells usi...



Western blot analysis of HepG2 cells usi...



Immunofluorescence analysis of HeLa cell...