

Product Datasheet

NFkB p65 (Phospho-Ser536) antibody (orb223069)

Description	Rabbit polyclonal antibody to NFkB p65 (Phospho-Ser536)
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	C-terminal
Concentration	1mg/ml
Note	For research use only
Clonality	Polyclonal
Uniprot ID	Q04206
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:40000
Expiration Date	12 months from date of receipt.