



Product Datasheet

SP1 (Phospho-Thr739) antibody (orb223101)

Description	Rabbit polyclonal antibody to SP1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P08047
Dilution Range	WB: 1:500-1:1000, ELISA: 1:5000
Expiration Date	12 months from date of receipt.



Western blot analysis of 3T3 cells using...



Western blot analysis of COS cells using...