



Product Datasheet

STAT3 (Phospho-Tyr705) antibody (orb223107)

Description	Rabbit polyclonal antibody to STAT3
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IHC, IP, WB
Immunogen	C-terminal
Concentration	1mg/ml
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Clonality	Polyclonal
Uniprot ID	P40763
Dilution Range	WB: 1:500-1:3000, IHC-P: 1:50-1:100, ELISA: 1:20000
Expiration Date	12 months from date of receipt.



Western blot analysis of HeLa cells using STAT-3 (pTyr705) antibody.



Immunohistochemical analysis of human brain tissue showing STAT3 (pTyr705) staining.