

---

## Product Datasheet

### Rabbit Human Ig Kappa Light Chain (free) antibody (orb22351)

**Description**

Precipitating polyclonal rabbit

chain. To detect and identify free ...

**Species/Host**

Rabbit

**Conjugation**

Unconjugated

**Preservatives**

Delipidated, heat inactivated lyophilized whole serum. No preservative added, as it may interfere with the antibody activity.

**Form/Appearance**

Delipidated, heat inactivated lyophilized whole serum. No preservative added, as it may interfere with the antibody activity.

**Storage**

Storage: The lyophilized product is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

A pool of purified Bence Jones kappa proteins isolated from human urine. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUbio accepts no liability for any inaccuracies or omissions in this information.

**Expiration Date**

12 months from date of receipt.