

Product Datasheet

6X His tag antibody (HRP) (orb223636)

Description	Mouse monoclonal antibody against 6X His tag conjugated to HRP
Species/Host	Mouse
Conjugation	HRP
Tested Applications	ELISA, IP, WB
Immunogen	PAX6 transcription factor linked to histidine tag
Target	6X His tag
Preservatives	Purified IgG conjugated to Horseradish Peroxidase using Innova Biosciences Lightning-Link, containing Phosphate buffered saline and 0.01% Thiomersal
Concentration	1mg/ml
Note	For research use only
Application notes	Cells from immunized Balb/c mice were fused with the NS1 myeloma cell line
Isotype	IgG1
Clonality	Monoclonal
Dilution Range	Optimal antibody dilution should be determined by titration however as a guideline use at 0.5ug-2ug/ml for WB
Expiration Date	12 months from date of receipt.