

---

## Product Datasheet

### Rabbit Mouse IgG Fab antibody (orb22377)

**Description**

Precipitating polyclonal rabbit antiserum

precipitating techniques a...

**Species/Host**

Rabbit

**Conjugation**

Unconjugated

**Preservatives**

Delipidated, heat inactivated lyophilized stable whole serum No preservative added, as it may interfere with the antibody activity. Total protein and IgG concentration in the antiserum are comparable to those of pooled rabbit serum. No foreign proteins ad

**Form/Appearance**

Delipidated, heat inactivated lyophilized stable whole serum No preservative added, as it may interfere with the antibody activity. Total protein and IgG concentration in the antiserum are comparable to those of pooled rabbit serum. No foreign proteins ad

**Storage**

Storage: The lyophilized antiserum is shipped at ambient temperature and may be stored at +4°C; prolonged storage at or below -20°C. Working dilutions are prepared by adding sterile phosphate buffered saline (PBS, pH 7.2). Repeated thawing and freezing should be avoided. Working dilutions should be stored at +4°C, not refrozen, and preferably used the same day. If a slight precipitation occurs upon storage, this should be removed by centrifugation. It will not affect the performance of the product

**Note**

For research use only

**Clonality**

Polyclonal

**Source**

Purified Fab of normal mouse IgG. Freund's complete adjuvant is used in the first step of the immunization procedure.

**Hazard Information**

This product is intended FOR RESEARCH USE ONLY, and FOR TESTS IN VITRO, not for use in diagnostic or therapeutic procedures involving humans or animals. This datasheet is as accurate as reasonably achievable, but Nordic-MUBio accepts no liability for any inaccuracies or omissions in this information.