

Product Datasheet

EpCAM antibody (FITC) (orb223770)

Description	Rat monoclonal antibody against EpCAM conjugated to FITC
Species/Host	Rat
Reactivity	Mouse
Conjugation	FITC
Tested Applications	ICC, IHC-Fr
Immunogen	Thymic epithelial cell line
Target	EpCAM
Preservatives	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
Concentration	1mg/ml
Note	For research use only
Application notes	X63Ag8.653
Isotype	IgG2a
Clonality	Monoclonal
Uniprot ID	Q99JW5
Entrez	17075
Dilution Range	Optimal antibody dilution should be determined by titration.
Expiration Date	12 months from date of receipt.