


## Product Datasheet

### F4/80 antibody (FITC) (orb223773)

<b>Description</b>	Gallus monoclonal antibody against F4/80 conjugated to FITC
<b>Species/Host</b>	Gallus
<b>Reactivity</b>	Mouse
<b>Conjugation</b>	FITC
<b>Tested Applications</b>	ICC
<b>Immunogen</b>	A synthetic peptide derived from C-terminus of Mouse F4/80
<b>Target</b>	F4/80
<b>Preservatives</b>	Purified IgG conjugated to Fluorescein using Innova Biosciences Lightning-Link, supplied in Phosphate buffered saline (PBS) containing 0.09% Sodium azide
<b>Concentration</b>	1mg/ml
<b>Note</b>	For research use only
<b>Isotype</b>	IgY
<b>Clonality</b>	Polyclonal
<b>Uniprot ID</b>	Q3S4B0
<b>Entrez</b>	<b>13733</b>
<b>Dilution Range</b>	Optimal antibody dilution should be determined by titration
<b>Expiration Date</b>	12 months from date of receipt.



Immunofluorescence analysis of using F4/...