

Product Datasheet

Human IL-19 protein (orb223823)

Description	Recombinant Human IL-19 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human IL-19
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder. Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and HPLC analyses
Source	E. coli
MW	17.9 kDa
Expiration Date	6 months from date of receipt.