


Product Datasheet

Human IL-36 beta protein (orb223824)

Description	Recombinant Human IL-36 beta protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human IL-36 beta
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder. Lyophilized from a 0.2 µm filtered concentrated solution in 2 x PBS, pH 7.4.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses
Source	E. coli
MW	17.2 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human IL-36 beta pr...