


## Product Datasheet

### Human sRANK Receptor protein (orb223829)

<b>Description</b>	Recombinant Human sRANK Receptor protein
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	SDS-PAGE
<b>Target</b>	Human sRANK Receptor
<b>Preservatives</b>	Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.0, 150mM NaCl.
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	> 98% by SDS-PAGE and HPLC analyses
<b>Source</b>	E. coli
<b>MW</b>	19.1 kDa
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human sRANK Recepto...