

Product Datasheet

Human HVEM-Fc protein (orb223830)

Description	Recombinant Human HVEM-Fc protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human HVEM-Fc
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2 µm filtered concentrated solution in PBS, pH7.4, with 3% Trehalose
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 95% by SDS-PAGE and HPLC analyses
Source	Pichia. Pastoris
MW	58 kDa
Expiration Date	6 months from date of receipt.