

Product Datasheet

Human NOV protein (orb223838)

| | |
|----------------------------|---|
| Description | Recombinant Human NOV protein |
| Reactivity | Human |
| Conjugation | Unconjugated |
| Tested Applications | SDS-PAGE |
| Target | Human NOV |
| Preservatives | Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.6, 150 mM NaCl |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only |
| Purity | > 95% by SDS-PAGE and HPLC analyses |
| MW | 36.2 kDa |
| Expiration Date | 6 months from date of receipt. |