

Product Datasheet

Human IGF-2 protein (orb223839)

Description	Recombinant Human IGF-2 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human IGF-2
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2 µm filtered concentrated solution in 20 mM Tris-HCl, pH 8.0, 150 mM NaCl, with 0.02% Tween-20.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 98% by SDS-PAGE and HPLC analyses
MW	7.5 kDa
Expiration Date	6 months from date of receipt.