


Product Datasheet

Human NRG1-beta 1 protein (orb223843)

Description	Recombinant Human NRG1-beta 1 protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human NRG1-beta 1
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder, Lyophilized from a 0.2µm filtered solution in 1 x PBS, pH 7.4, with 5% trehalose
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses
MW	7.5 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human NRG1-beta 1 p...