

## Product Datasheet

### Rat HB-EGF protein (orb223887)

<b>Description</b>	Rat HB-EGF protein
<b>Reactivity</b>	Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	SDS-PAGE
<b>Target</b>	Rat HB-EGF
<b>Preservatives</b>	Sterile Filtered White lyophilized (freeze-dried) powder.
<b>Storage</b>	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	> 95% by SDS-PAGE and HPLC analyses.
<b>Source</b>	E. coli
<b>MW</b>	9.7 kDa
<b>Expiration Date</b>	6 months from date of receipt.