


Product Datasheet

Equine IL-2 Cys protein (orb223896)

Description	Recombinant Equine IL-2 Cys protein
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Equine IL-2 Cys
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 97% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	14.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Equine IL-2 Cys pro...