


Product Datasheet

Murine MIP-1 gamma protein (orb223915)

| | |
|----------------------------|--|
| Description | Recombinant Murine MIP-1 gamma protein |
| Conjugation | Unconjugated |
| Tested Applications | SDS-PAGE |
| Target | Murine MIP-1 gamma |
| Preservatives | Sterile Filtered White lyophilized (freeze-dried) powder. |
| Storage | Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles. |
| Note | For research use only |
| Purity | > 95% by SDS-PAGE and HPLC analyses. |
| Source | E. coli |
| MW | 11.6 kDa |
| Expiration Date | 6 months from date of receipt. |



SDS-PAGE analysis of Murine MIP-1 gamma ...