

Product Datasheet

Lamin A/C Antibody (Phospho-Ser392) (orb2239447)

Description	Lamin A/C Antibody (Phospho-Ser392)
Species/Host	Rabbit
Conjugation	Unconjugated
Tested Applications	ELISA, ICC, IF, IHC, WB
Immunogen	A synthesized peptide derived from human Lamin A/C around the phosphorylation site of Serine 392
Preservatives	Liquid PBS (pH 7.4) 150mM NaCl, 0.02% sodium azide and 50% glycerol
Form/Appearance	Liquid PBS (pH 7.4) 150mM NaCl, 0.02% sodium azide and 50% glycerol
Storage	Store at -20°C. Stable for 12 months from date of receipt.
Note	For research use only
Application notes	Application Info: WB: 1:500~2000IHC: 1:50~1000IF/ICC: 1:100~500P-ELISA: 1:20000~40000
Isotype	IgG
Clonality	Polyclonal
Source	Rabbit
MW	74 kDa, 65 kDa
Uniprot ID	P02545
Expiration Date	12 months from date of receipt.