


Product Datasheet

Human DHH C23II protein (orb223945)

Description	Recombinant Human DHH C23II protein
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	SDS-PAGE
Target	Human DHH C23II
Preservatives	Sterile Filtered White lyophilized (freeze-dried) powder.
Storage	Store at 4°C for up to two weeks. For long term storage, aliquot and store at -20°C, avoid freeze/thaw cycles.
Note	For research use only
Purity	> 96% by SDS-PAGE and HPLC analyses.
Source	E. coli
MW	19.9 kDa
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human DHH C23II pro...