

Product Datasheet

MKK3/6 Peptide (orb224188)

Description	The peptide is used to block MKK3/6 antibody (orb223989) reactivity
Reactivity	Human, Mouse, Porcine, Rat, Zebrafish
Tested Applications	Blocking
Target	MAP2K3; MAP2K6
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb223989. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P46734 , O09110
Entrez	5606 , 26397
Expiration Date	12 months from date of receipt.