

Product Datasheet

S6K1 (phospho-S447) Peptide (orb224201)

Description	The peptide is used to block S6K1 (phospho-S447) antibody (orb224002) reactivity
Reactivity	Human, Mouse, Porcine, Rabbit, Rat
Tested Applications	Blocking
Target	RPS6KB1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224002. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P67999 , P23443 , Q8BSK8
Entrez	6198 , 72508 , 83840
Expiration Date	12 months from date of receipt.