

Product Datasheet

STAT4 (phospho-Y693) Peptide (orb224210)

Description	The peptide is used to block STAT4 (phospho-Y693) antibody (orb224011) reactivity
Reactivity	Human, Mouse, Porcine, Rabbit, Rat, Sheep
Tested Applications	Blocking
Target	STAT4
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224011. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P42228 , Q14765
Entrez	6775 , 20849
Expiration Date	12 months from date of receipt.