

Product Datasheet

Synapsin 1 (phospho-S9) Peptide (orb224218)

Description	The peptide is used to block Synapsin 1 (phospho-S9) antibody (orb224019) reactivity
Reactivity	Human, Mouse, Porcine, Rat
Tested Applications	Blocking
Target	SYN1
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224019. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	O88935 , P17600 , P09951
Entrez	20964 , 6853 , 24949
Expiration Date	12 months from date of receipt.