

Product Datasheet

Histone H4 (phospho-S47) Peptide (orb224229)

Description	The peptide is used to block Histone H4 (phospho-S47) antibody (orb224030) reactivity
Reactivity	Human, Rat
Tested Applications	Blocking
Target	HIST1H4A
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224030. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P62805
Entrez	8368, 8359, 8365, 8360, 8364, 8366, 121504, 8363, 8367, 8362, 554313, 8361, 8370, 8294
Expiration Date	12 months from date of receipt.