

Product Datasheet

TSC2 (phospho-T1462) Peptide (orb224237)

Description	The peptide is used to block TSC2 (phospho-T1462) antibody (orb224038) reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	TSC2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224038. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P49816 , P49815 , Q61037
Entrez	7249 , 24855
Expiration Date	12 months from date of receipt.