

Product Datasheet

MOV10L1 Peptide (orb224303)

| | |
|----------------------------|--|
| Description | The peptide is used to block MOV10L1 antibody (orb224104) reactivity |
| Reactivity | Human, Mouse, Rat |
| Tested Applications | Blocking |
| Target | MOV10L1 |
| Form/Appearance | Lyophilized powder |
| Storage | Shipped at 4°C. Store at -20°C for one year. |
| Note | For research use only |
| Application notes | This peptide was used as immunogen to produce orb224104. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours. |
| Purity | > 85% |
| Source | Synthetic |
| Uniprot ID | Q9BXT6 , Q99MV5 |
| Entrez | 83456 , 54456 |
| Expiration Date | 12 months from date of receipt. |