

Product Datasheet

Histone Deacetylase 8 (phospho-S39) Peptide (orb224306)

Description	The peptide is used to block Histone Deacetylase 8 (phospho-S39) antibody (orb224107) reactivity
Reactivity	Human, Mouse, Primate, Rat
Tested Applications	Blocking
Target	HDAC8
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224107. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q9BY41 , B1WC68 , Q8VH37
Entrez	70315 , 363481 , 100911968 , 55869
Expiration Date	12 months from date of receipt.