

Product Datasheet

EEF2K (phospho-S366) Peptide (orb224354)

Description	The peptide is used to block EEF2K (phospho-S366) antibody (orb224155) reactivity
Reactivity	Human, Primate, Rat
Tested Applications	Blocking
Target	EEF2K
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224155. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	P70531 , O00418
Entrez	101930123 , 29904 , 25435
Expiration Date	12 months from date of receipt.