

## Product Datasheet

# Alpha-adducin (phospho-S726/713) Peptide (orb224356)

<b>Description</b>	The peptide is used to block Alpha-adducin (phospho-S726/713) antibody (orb224157) reactivity
<b>Reactivity</b>	Human, Mouse, Porcine, Rat
<b>Tested Applications</b>	Blocking
<b>Target</b>	ADD1; ADD2
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb224157. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">Q63028</a> , <a href="#">P35612</a> , <a href="#">P35611</a> , <a href="#">Q9QYB8</a> , <a href="#">Q05764</a> , <a href="#">Q9QYCO</a>
<b>Entrez</b>	<a href="#">11519</a> , <a href="#">119</a> , <a href="#">118</a> , <a href="#">24170</a> , <a href="#">24171</a> , <a href="#">11518</a>
<b>Expiration Date</b>	12 months from date of receipt.