

## Product Datasheet

### Ah Receptor (phospho-S36) Peptide (orb224363)

<b>Description</b>	The peptide is used to block Ah Receptor (phospho-S36) antibody (orb224164) reactivity
<b>Reactivity</b>	Human, Mouse, Rabbit, Rat, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	AHR
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb224164. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P30561</a> , <a href="#">A9YTQ3</a> , <a href="#">P41738</a> , <a href="#">Q3U1U7</a> , <a href="#">Q75NT5</a> , <a href="#">P35869</a>
<b>Entrez</b>	<a href="#">11624</a> , <a href="#">57491</a> , <a href="#">196</a> , <a href="#">25690</a> , <a href="#">11622</a> , <a href="#">498999</a>
<b>Expiration Date</b>	12 months from date of receipt.