

Product Datasheet

FOXC1/2 Peptide (orb224368)

Description	The peptide is used to block FOXC1/2 antibody (orb224169) reactivity
Reactivity	Human, Mouse, Rat
Tested Applications	Blocking
Target	FOXC1; FOXC2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224169. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q99958 , Q12948 , Q63246 , Q61572 , Q61850
Entrez	2303 , 14234 , 17300 , 2296
Expiration Date	12 months from date of receipt.