

## Product Datasheet

### TFRC Antibody (Phospho-Ser24) (orb2243687)

<b>Description</b>	TFRC Antibody (Phospho-Ser24)
<b>Species/Host</b>	Rabbit
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ELISA, IF, IHC, IHC-P, WB
<b>Immunogen</b>	The antiserum was produced against synthesized peptide derived from human CD71/TfR around the phosphorylation site of Ser24.
<b>Preservatives</b>	Liquid
<b>Form/Appearance</b>	Liquid
<b>Storage</b>	-20°C
<b>Note</b>	For research use only
<b>Application notes</b>	Application Info: WB: 1:500~1:3000; IHC: 1:50~1:100; ELISA: 1:10000. Dilution Info: The optimal working dilution should be determined by the end user.
<b>Clonality</b>	Polyclonal
<b>MW</b>	84 kDa
<b>Uniprot ID</b>	<a href="#">P02786</a>
<b>NCBI</b>	<a href="#">NP_001121620.1</a>
<b>Expiration Date</b>	12 months from date of receipt.