

## Product Datasheet

# IGF1 Receptor (phospho-Y1161) Peptide (orb224370)

<b>Description</b>	The peptide is used to block IGF1 Receptor (phospho-Y1161) antibody (orb224171) reactivity
<b>Reactivity</b>	Human, Mouse, Rat, Zebrafish
<b>Tested Applications</b>	Blocking
<b>Target</b>	IGF1R
<b>Form/Appearance</b>	Lyophilized powder
<b>Storage</b>	Shipped at 4°C. Store at -20°C for one year.
<b>Note</b>	For research use only
<b>Application notes</b>	This peptide was used as immunogen to produce orb224171. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
<b>Purity</b>	> 85%
<b>Source</b>	Synthetic
<b>Uniprot ID</b>	<a href="#">P24062</a> , <a href="#">Q60751</a> , <a href="#">P15208</a> , <a href="#">P08069</a> , <a href="#">P15127</a> , <a href="#">P06213</a>
<b>Entrez</b>	<a href="#">3480</a> , <a href="#">16337</a> , <a href="#">16001</a> , <a href="#">25718</a> , <a href="#">3643</a>
<b>Expiration Date</b>	12 months from date of receipt.