

Product Datasheet

MKK1/2 (phospho-S218/222) Peptide (orb224375)

Description	The peptide is used to block MKK1/2 (phospho-S218/222) antibody (orb224176) reactivity
Reactivity	Human, Mouse, Rabbit, Rat, Zebrafish
Tested Applications	Blocking
Target	MAP2K1; MAP2K2
Form/Appearance	Lyophilized powder
Storage	Shipped at 4°C. Store at -20°C for one year.
Note	For research use only
Application notes	This peptide was used as immunogen to produce orb224176. Experiment Notes: Blocking Peptide to the diluted primary antibody in a molar ratio of 10:1 (peptide to antibody) and incubate the mixture at 4°C for overnight or at room temperature for 2 hours.
Purity	> 85%
Source	Synthetic
Uniprot ID	Q02750 , Q01986 , P36506 , P31938 , Q63932 , P36507
Entrez	26395 , 170851 , 58960 , 5605 , 5604 , 26396
Expiration Date	12 months from date of receipt.